SISCO-4 in 1 Digital Water Quality Tester

sisco

Please read this manual carefully before use

Technical Parameters

Operation Manual
C:0 M:66 Y

195mm

model	SISCO-DDS-10		
Conductivity	Measuring range	0.0µS/cm(ppm)-20.00mS/cm(ppt)	
	resolution	0.1µS/cm(ppm)-0.01mS/cm(ppt)	
	Measuring accuracy	±1%F.S	
	Measuring range	0.0-60ÿ/32-140ÿ	
	resolution	0.1℃	
	Compensation	Automatic, 0-60 C/32-140 F	
temperature	range Temperature	10-40°C (default: 25°C)	
	coefficient Reference	± 1 % F. S	
	temperature range	0-20.00ppt	
	Measuring accuracy Temperature	Extraposition	
TDS	sensor structure Measuring range	External 0-20.00ppt	
	TDS coefficient	0.01-1.0 adjustable (default 0.55)	
	measurement	% F. S	
	accuracy	0-16.00ppt NACL	
salinity	measurement	NACL 0.65	
	range salinity	± 1 % F. S	
	coefficient measurement	U-shaped glass	
other	accuracy electrode	Auto; 1408µS/cm (1 point calibration)	
	sensor type	rice digital screen	
	calibration mode	orange	
	Sound switch	Possess	
	backlight color sound switching pow	(7#)AAA battery X1PCS	

Troubleshooting

error code	Fault diagnosis	Solution
E-01	The electrode sensor is damaged or aged	Replace electrode tip
E-02	Battery power is too low	Replace the battery, it is recommended to use (7#) AAA alkaline batteries
E-03	Unknown fault	Free replacement for non-human reasons during the warranty period
E-04	beyond the scope of use	Clean the electrode sensor and use it according to the instructions
E-05	Other hardware failures	Free replacement for non-human reasons during the warranty period

- 1 -

Y:91 K:0

130mm

packing list


Host×1

Calibration solution ×1

Instructions, certificate ×1

(AAA)7# alkaline battery×1

product structure



sisco

Key Function

button	Function
\bigcirc	 Long press to turn on and off, short press to turn on or off the backlight.
	 Press and hold to enter the calibration mode setting state, select the current option to enter step adjustment
(SET	 Enter the setting state. In the setting state, cycle into the next setting menu setting state. Adjust the current menu option. In the setting state, adjust the step.
ENT	 In the calibration mode, save the calibration and exit the calibration mode. In the setting state, save the current setting options or in the step setting state, exit the setting menu.

Settings menu description

Show menu	Menu function description	Settable options Opt	ion function description Fa	ctory default
UNIT	temperature scale unit	ÿ	Celsius	ÿ
		ÿ	Fahrenheit	
LOCK	Auto lock	O N	On Off	ÿ
		OFF	On Off	
BUZZ	Buzzer switch	O N	Step	ÿ
		OFF	0.01	
ТС	Temperature coefficient	0.00-0.05	Step 1 Step	0.02
MRT	reference temperature	10-40	0.01	2 5
TDSc	TDS coefficient	0.01-1.00		0.55
HALL	salinity coefficient	NACL 0.65		ÿ
rESE	reset	No	No	
		Yes	restart	

- 3 -

sisco

Conductivity meter calibration

Calibration is required when using it for the first time, replacing the battery, or replacing the electrode head. Regular calibration is

recommended. The standard solution should be prepared and used as soon as possible. Do not reuse it. Standard solutions stored improperly or for too long will affect the measurement accuracy.

Rinse in advance to improve calibration accuracy.

 Prepare 1408yS/cm standard solution. Press and hold the power button to turn on the meter and make sure the conductivity meter is in the conductivity measurement mode (the screen displays ÿS/cm), otherwise press the button to switch to the Hold under the conductivity measurement mode. 	uring and the second s
3. Clean the electrode sensor with pure water, and completely immerse the electrode sensor in the 1408 ÿS/cm standard solution.	
Press and hold the calibration key, and the screen measurement area will display and flash, which measurement area will display	
the calibration process. At this time, 1408 will be displayed on the lower right corner of the screen, gently stir the	17855
conductivity meter a few times, wait for the screen to display PASS, and press the confirmation key.	µS/cm Offset
Press the ENT key, calibration is completed.	1408

Instructions

195mm

Pre-rinsing the electrode sensor every time it is placed into the sample to be measured can improve test accuracy.





mode



1. Press and hold the power button to un on the

device, and remove the conductivity meter electrode cap. cm for c

2. Press the bar on to switch to the desired measurement mode. The screen parameter display area displays ÿS/

cm for conductivity measurement mode, ppm for TDS measurement mode, and ppt for salinity measurement





3. Completely immerse the electrode sensor into the water sample to be measured, and stir it gently a few times. After the measured value stabilizes icon is displayed at the top of the

screen, the current val. Sisplayed on the screen is the measured value of the sample to be measured, and the value of the water sample to be measured is displayed in the

lower right corner. temperature.

4. After use, please clean the electrode sensor with pure water.

Note: When measuring, please be sure to confirm whether the current is the required test mode!

Replacement battery

1. Turn left to remove the top cover

2. Push open the battery compartment

cover 3. Take out the old battery and put in the new battery

(It is recommended to use 7# AAA alkaline dry batteries)

4. Close the battery compartment

cover 5. Turn right to tighten the top cover



Replace electrode head

1. Pull off the electrode cap directly

2. Unplug the electrode head, do not rotate it to avoid damaging the circuit. 3. Align the two electrode

slots and push the electrode head into the main unit. 4. Tighten the electrode fixing ring by turning

right. 5. Cover the electrode cap.



- 4 -

- 5 -

Note: During installation, please confirm that the waterproof ring in the electrode fixing ring is reset, and please tighten the electrode fixing ring with appropriate force to avoid affecting the waterproof and dustoroof effect.

Note: When installing, please pay attention to whether the waterproof ring in the top cover is reset, and please tighten the top cover with appropriate force to avoid affecting the waterproof and dustproof effect.